

Product Information

MemDX™ Membrane Protein Human XKR7 (XK related 7) Expressed *in vitro* *E.coli* expression system, Full Length

Cat. No.: **MPX3640K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human XKR7 membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

XKR7

Protein Length

Full Length

Protein Class

Receptor

TMD

7

Sequence

MAAKSDGAAASASPDPEGAAGGARGSAGGRGEAAAAAGPPGVVAGGGPGPRYELRDCCWVLCALLVFFSDGATDLWLAASYL

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

XKR7

Full Name

XK related 7

Alternative Names

XKR7; C20orf159; dJ310O13.4; XK-related protein 7; X Kell blood group precursor-related family, member 7; X-linked Kx blood group related 7; XK, Kell blood group complex subunit-related family, member 7; XK related 7

Gene ID

[343702](#)

UniProt ID

[Q5GH72](#)