

# Product Information

## **MemDX™ Membrane Protein Mouse Snn (Stannin) Expressed *in vitro* E.coli expression system, Full Length**

Cat. No.: **MPX3674K**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse Snn membrane protein expressed *in vitro* E.coli expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### **Host Species**

Mouse

#### **Target Protein**

Snn

#### **Protein Length**

Full Length

#### **Protein Class**

Receptor

#### **TMD**

1

#### **Sequence**

MSIMDHSPPTTGVVTVIVILIAIALGALILGCWCYLRLQRISQSEDEESIVGDGETKEPFLLVQYSAKGPCVERKAKLMTANSPEVHG

### Product Description

#### **Expression Systems**

*in vitro* E.coli expression system

#### **Tag**

10xHis tag at the N-terminus

#### **Protein Format**

Soluble

#### **Form**

Liquid or Lyophilized powder

#### **Buffer**

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

Snn

**Full Name**

Stannin

**Alternative Names**

Snn; AI848521; AW547286; 2810407J07Rik; Stannin

**Gene ID**

[20621](#)

**UniProt ID**

[P61807](#)