

Product Information

MemDX™ Membrane Protein Human MFAP3L (Microfibril associated protein 3 like)

Expressed *in vitro* *E.coli* expression system, Full Length of Mature Protein

Cat. No.: **MPX3871K**

This product is for research use only and is not intended for diagnostic use.

This product is a Human MFAP3L membrane protein expressed *in vitro* *E.coli* expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MFAP3L

Protein Length

Full Length of Mature Protein

Protein Class

Receptor

TMD

1

Sequence

KSVTNSTLNGTNVVLGSVPVVIARTDHIIVKEGNSALINCSVYGIPDPQFKWYNSIGKLLKEEEDKERGGGKWQMHD SGLLNITKVS

Product Description

Expression Systems

in vitro *E.coli* expression system

Tag

10xHis tag at the N-terminus

Protein Format

Soluble

Form

Liquid or Lyophilized powder

Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

MFAP3L

Full Name

Microfibril associated protein 3 like

Alternative Names

MFAP3L; NYD-sp9; microfi brillar-associated protein 3-like; microfibrillar associated protein 3 like; testis development protein NYD-SP9; Microfibril associated protein 3 like

Gene ID

[9848](#)

UniProt ID

[O75121](#)