

Product Information

MemDX™ Human VTCN1 HEK293T-OS8 Cell Line

Cat. No.: **S01YF-0123-KX267**

This product is for research use only and is not intended for diagnostic use.

Product Information

Target Protein

VTCN1

Target Protein Species

Human

Accession Number

NM_024626.3

Protein Tag

Tag-free

Host Cell Type

HEK293T-OS8

Target Classification

Immune Checkpoint

Target Family

Immune Checkpoint

Target Research Area

Immunology Research

Product Properties

Morphology

Fibroblastoid cells growing as a monolayer

Assay Types

Drug screening and biological assays

Resistance

Puromycin; Hygromycin

Mycoplasma Testing

Negative

Sterility Testing

10 passages

Biosafety Level

Level 1

Activity

Yes

Quantity

5x10⁶ cells

Form

Frozen cells

Freeze Medium

70% DMEM + 20% FBS + 10% DMSO

Culture Medium

DMEM + 10% FBS + 0.5µg/mL Puromycin + 50µg/mL Hygromycin

Selective Antibiotic(s)

Regular antibiotics active against mycoplasmas, bacteria and fungi.

Handling Notes

Frozen cells should be thawed immediately upon receipt and grown according to handling procedure to ensure cell viability and proper assay performance.

Note: Do not freeze the cells upon receipt as it may result in irreversible damage to the cell line.

Disclaimer: We cannot guarantee cell viability if the cells are not thawed immediately upon receipt and grown according to handling procedure.

Incubation

37°C with 5% CO₂

Applications

Drug screening and biological assays

Application Notes

Cells were plated in a 384-well plate and incubated overnight at 37°C and 5% CO₂ to allow the cells to attach and grow. Cells were then stimulated with a control for high-throughput drugs screening and functional assays.

Use Restrictions

These cells are distributed for research use only.

Shipping

Dry ice

Storage

Liquid nitrogen

Target**Full Name**

V-set domain containing T cell activation inhibitor 1

Introduction

This gene encodes a protein belonging to the B7 costimulatory protein family. Proteins in this family are present on the surface of antigen-presenting cells and interact with ligand bound to receptors on the surface of T cells. Studies have shown that high levels of the encoded protein has been correlated with tumor progression. A pseudogene of this gene is located on chromosome 20. Multiple transcript variants encoding different isoforms have been found for this gene.

Alternative Names

VTCN1; B7X; B7H4; B7S1; B7-H4; B7h.5; VCTN1; PRO1291; V-set domain-containing T-cell activation inhibitor 1; B7 family member, H4; B7 homolog 4; B7 superfamily member 1; T cell costimulatory molecule B7x; immune costimulatory protein B7-H4; V-set domain containing T cell activation inhibitor 1

Gene ID

[79679](#)

UniProt ID

[Q7Z7D3](#)