

Product Information

mPEG-DSPE, 2000

Phosphatidylethanolamine

Cat. No.: **S01YF-0324-KX398**

This product is for research use only and is not intended for diagnostic use.

It is a phospholipid related product for research use only.

Product Type

Phospholipids

Classification

Phosphatidylethanolamine (PE)

Appearance

Solid

CAS Number

474922-77-5

Molecular Weight

~2820 (mPEG MW = 2000 Da)

Purity

>95%

Size

50mg; 100mg; 250mg

Storage

-20°C or below