

Product Information

MemDX™ Human TACSTD2 CT26 Cell Line

Cat. No.: **S01YF-0324-KX67**

This product is for research use only and is not intended for diagnostic use.

Product Information

Target Protein

TACSTD2

Target Protein Species

Human

Host Cell Type

CT26

Target Classification

Others

Target Family

Others

Target Research Area

Ocular Research

Related Diseases

Corneal Dystrophy

Product Properties

Assay Types

Functional assay and biological assay

Mycoplasma Testing

Negative

Biosafety Level

Level 1

Activity

Yes

Form

Frozen cells

Freeze Medium

DMEM + 20% FBS + 10% DMSO

Culture Medium

DMEM + 10% FBS

Selective Antibiotic(s)

Regular antibiotics active against mycoplasmas, bacteria and fungi.

Handling Notes

Frozen cells should be thawed immediately upon receipt and grown according to handling procedure to ensure cell viability and proper assay performance.

Note: Do not freeze the cells upon receipt as it may result in irreversible damage to the cell line.

Disclaimer: We cannot guarantee cell viability if the cells are not thawed immediately upon receipt and grown according to handling procedure.

Incubation

37°C with 5% CO₂

Applications

Drug screening and biological assays

Application Notes

Cells were plated in a 384-well plate and incubated overnight at 37°C and 5% CO₂ to allow the cells to attach and grow. Cells were then stimulated with a control for high-throughput drugs screening and functional assays.

Use Restrictions

These cells are distributed for research use only.

Shipping

Dry ice

Storage

Liquid nitrogen

Target

Full Name

Tumor associated calcium signal transducer 2

Introduction

This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

Alternative Names

TACSTD2; EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1; 40kD glycoprotein, identified by monoclonal antibody GA733; cell surface glycoprotein TROP2; cell surface glycoprotein Trop-2; epithelial glycoprotein-1; gastrointestinal tumor-associated antigen GA7331; membrane component, chromosome 1, surface marker 1; pancreatic carcinoma marker protein GA733-1; pancreatic carcinoma marker protein GA7331; trophoblast cell surface antigen 2; truncated TACSTD2; Tumor associated calcium signal transducer 2

Gene ID

4070

UniProt ID

P09758