

Product Information

MemDX™ Membrane Protein Human SIGLEC10 (17-546aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX108**

This product is for research use only and is not intended for diagnostic use.

This product is a Human SIGLEC10 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

SIGLEC10

Protein Length

Partial (17-546aa)

Molecular Weight

59.0 kDa

Sequence

17-546aa (Accession # Q96LC7)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

SIGLEC10

Full Name

Sialic acid binding Ig like lectin 10

Alternative Names

SIGLEC10; SLG2; PRO940; SIGLEC-10; sialic acid-binding Ig-like lectin 10; sialic acid binding Ig-like lectin 10 Ig-like lectin 7; siglec-like gene 2; siglec-like protein 2; Sialic acid binding Ig like lectin 10

Gene ID

[89790](#)

UniProt ID

[Q96LC7](#)