

Product Information

MemDX™ Membrane Protein Human CCR4 (1-39; 99-111; 176-200; 268-284aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX126**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CCR4 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CCR4

Protein Length

Partial (1-39; 99-111; 176-200; 268-284aa)

Molecular Weight

37.1 kDa

Sequence

1-39; 99-111; 176-200; 268-284aa (Accession # P51679)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CCR4

Full Name

C-C motif chemokine receptor 4

Alternative Names

CKR4; K5-5; CD194; CMKBR4; ChemR13; CC-CKR-4; HGCN:14099; C-C chemokine receptor type 4; C-C CKR-4; CCR-4; chemokine (C-C motif) receptor 4; chemokine (C-C) receptor 4; CCR4; C-C motif chemokine receptor 4

Gene ID

[1233](#)

UniProt ID

[Q9ULM6](#)