

Product Information

MemDX™ Membrane Protein Human GPC3 (511-560aa) Expressed in HEK293, N-Human Fc tag

Cat. No.: **S01YF-0423-KX127**

This product is for research use only and is not intended for diagnostic use.

This product is a Human GPC3 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

GPC3

Protein Length

Partial (511-560aa)

Molecular Weight

43.3 kDa

Sequence

511-560aa (Accession # P51654)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

N-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

GPC3

Full Name

Glypican 3

Alternative Names

SGB; DGSX; MXR7; SDYS; SGBS; OCI-5; SGBS1; GTR2-2; glypican-3; glypican proteoglycan 3; heparan sulphate proteoglycan; intestinal protein OCI-5; secreted glypican-3

Gene ID

[2719](#)

UniProt ID

[P51654](#)