

Product Information

MemDX™ Membrane Protein Human IL2 (21-153aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX151**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL2 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL2

Protein Length

Partial (21-153aa)

Molecular Weight

41.6 kDa

Sequence

21-153aa (Accession # P60568)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

IL2

Full Name

Interleukin 2

Alternative Names

IL-2; TCGF; lymphokine; interleukin-2; T cell growth factor; aldesleukin; involved in regulation of T-cell clonal expansion

Gene ID

[3558](#)

UniProt ID

[P60568](#)