

Product Information

MemDX™ Membrane Protein Human TLR3 (28-80aa) Expressed in HEK293, C-Mouse Fc tag

Cat. No.: **S01YF-0423-KX162**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TLR3 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

TLR3

Protein Length

Partial (28-80aa)

Molecular Weight

103.5 kDa

Sequence

28-80aa (Accession # O15455)

Product Description

Expression Systems

HEK293

Tag

C-Mouse Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

TLR3

Full Name

Toll like receptor 3

Alternative Names

TLR3; CD283; IIAE2; Toll like receptor 3

Gene ID

[7098](#)

UniProt ID

[Q15455](#)