

# Product Information

## **MemDX™ Membrane Protein Human CD200 (31-232aa) Expressed in HEK293, C-6×His tag**

Cat. No.: **S01YF-0423-KX169**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CD200 membrane protein expressed in HEK293.

### Product Specifications

#### **Host Species**

Human

#### **Target Protein**

CD200

#### **Protein Length**

Partial (31-232aa)

#### **Molecular Weight**

23.3 kDa

#### **Sequence**

31-232aa (Accession # P41217)

### Product Description

#### **Activity**

Yes

#### **Expression Systems**

HEK293

#### **Tag**

C-6×His tag

#### **Protein Format**

Soluble

#### **Form**

LYOPH

#### **Purity**

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

#### **Buffer**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

**Storage**

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

**Target****Target Protein**

CD200

**Full Name**

CD200 molecule

**Alternative Names**

CD200; MRC; MOX1; MOX2; OX-2; OX-2 membrane glycoprotein; CD200 antigen; antigen identified by monoclonal antibody MRC OX-2; CD200 molecule

**Gene ID**

[4345](#)

**UniProt ID**

[P41217](#)