

Product Information

MemDX™ Membrane Protein Human CD200 (31-232aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX169**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CD200 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CD200

Protein Length

Partial (31-232aa)

Molecular Weight

23.3 kDa

Sequence

31-232aa (Accession # P41217)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CD200

Full Name

CD200 molecule

Alternative Names

CD200; MRC; MOX1; MOX2; OX-2; OX-2 membrane glycoprotein; CD200 antigen; antigen identified by monoclonal antibody MRC OX-2; CD200 molecule

Gene ID

[4345](#)

UniProt ID

[P41217](#)