

Product Information

MemDX™ Membrane Protein Human CD27 (20-191aa) Expressed in HEK293, C-Human Fc

Cat. No.: **S01YF-0423-KX176**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CD27 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CD27

Protein Length

Partial (20-191aa)

Molecular Weight

45.3 kDa

Sequence

20-191aa (Accession # P26842)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CD27

Full Name

CD27 molecule

Alternative Names

CD27; T14; S152; Tp55; TNFRSF7; S152. LPFS2; CD27 antigen; CD27L receptor; T cell activation antigen S152; T-cell activation antigen CD27; tumor necrosis factor receptor superfamily, member 7; CD27 molecule

Gene ID

[939](#)

UniProt ID

[P26842](#)