

Product Information

MemDX™ Membrane Protein SARS-CoV-2 S (310-560aa) Expressed in HEK293, C- 6×His tag

Cat. No.: **S01YF-0423-KX193**

This product is for research use only and is not intended for diagnostic use.

This product is a SARS-CoV-2 S membrane protein expressed in HEK293.

Product Specifications

Host Species

SARS-CoV-2

Target Protein

S

Protein Length

Partial (310-560aa)

Molecular Weight

31.2 kDa

Sequence

310-560aa (Accession # P0DTC2)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

S

Full Name

Surface glycoprotein

Alternative Names

S; structural protein; spike protein; Surface glycoprotein

Gene ID

[43740568](#)

UniProt ID

[P0DTC2](#)