

Product Information

MemDX™ Membrane Protein Mouse FAP (26-761aa) Expressed in HEK293, C-Mouse Fc tag

Cat. No.: **S01YF-0423-KX219**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse FAP membrane protein expressed in HEK293.

Product Specifications

Host Species

Mouse

Target Protein

FAP

Protein Length

Partial (26-761aa)

Molecular Weight

111.4 kDa

Sequence

26-761aa (Accession # P97321)

Product Description

Expression Systems

HEK293

Tag

C-Mouse Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

FAP

Full Name

Fibroblast activation protein

Alternative Names

SIMP; prolyl endopeptidase FAP; FAPalpha; dipeptidyl peptidase FAP; fibroblast activation protein alpha; gelatine degradation protease FAP; integral membrane serine protease; post-proline cleaving enzyme; seprase; serine integral membrane protease; surface-expressed protease

Gene ID

[14089](#)

UniProt ID

[P97321](#)