

Product Information

MemDX™ Membrane Protein Human IL6R (20-365aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX225**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL6R membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL6R

Protein Length

Partial (20-365aa)

Molecular Weight

64.7 kDa

Sequence

20-365aa (Accession # P08887)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

IL6R

Full Name

Interleukin 6 receptor

Alternative Names

IL6R; IL6Q; gp80; CD126; HIES5; IL-6R; IL6RA; IL6RQ; IL-1Ra; IL-6RA; IL6QTL; IL-6R-1; interleukin-6 receptor subunit alpha; CD126 antigen; IL-6 receptor subunit alpha; membrane glycoprotein 80; Interleukin 6 receptor

Gene ID

[3570](#)

UniProt ID

[P08887](#)