

Product Information

MemDX™ Membrane Protein Human NCAM1 (220-708aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX230**

This product is for research use only and is not intended for diagnostic use.

This product is a Human NCAM1 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

NCAM1

Protein Length

Partial (220-708aa)

Molecular Weight

102.3 kDa

Sequence

220-708aa (Accession # P13591)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 90% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

NCAM1

Full Name

Neural cell adhesion molecule 1

Alternative Names

NCAM1; CD56; NCAM; MSK39; antigen recognized by monoclonal antibody 5.1H11; neural cell adhesion molecule, NCAM; Neural cell adhesion molecule 1

Gene ID

[4684](#)

UniProt ID

[P13591](#)