

Product Information

MemDX™ Membrane Protein Human DLL3 (27-490aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX293**

This product is for research use only and is not intended for diagnostic use.

This product is a Human DLL3 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

DLL3

Protein Length

Partial (27-490aa)

Molecular Weight

74.4 kDa

Sequence

27-490aa (Accession # Q9NYJ7)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 80% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

DLL3

Full Name

Delta like canonical Notch ligand 3

Alternative Names

DLL3; SCDO1; delta-like protein 3; delta-like 3; delta3; drosophila Delta homolog 3; Delta like canonical Notch ligand 3

Gene ID

[10683](#)

UniProt ID

[Q9NYJ7](#)