

Product Information

MemDX™ Membrane Protein Human SDC1 (23-254aa) Expressed in HEK293, C-Human Fc and 6xHis tag

Cat. No.: **S01YF-0423-KX32**

This product is for research use only and is not intended for diagnostic use.

This product is a Human SDC1 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

SDC1

Protein Length

Partial (23-254aa)

Molecular Weight

50.8 kDa

Sequence

23-254aa (Accession # P18827)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc and 6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

SDC1

Full Name

Syndecan 1

Alternative Names

SDC1; SDC; CD138; SYND1; syndecan; syndecan-1; CD138 antigen; heparan sulfate proteoglycan fibroblast growth factor receptor; syndecan proteoglycan 1; Syndecan 1

Gene ID

[6382](#)

UniProt ID

[P18827](#)