

Product Information

MemDX™ Membrane Protein Human NOTCH3 (40-467aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX331**

This product is for research use only and is not intended for diagnostic use.

This product is a Human NOTCH3 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

NOTCH3

Protein Length

Partial (40-467aa)

Molecular Weight

70.9 kDa

Sequence

40-467aa (Accession # Q9UM47)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

NOTCH3

Full Name

Notch receptor 3

Alternative Names

NOTCH3; IMF2; LMNS; CASIL; CADASIL; CADASIL1; neurogenic locus notch homolog protein 3; Notch homolog 3; notch 3; Notch receptor 3

Gene ID

[4854](#)

UniProt ID

[Q9UM47](#)