

Product Information

MemDX™ Membrane Protein Human FZD10 (21-161aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX335**

This product is for research use only and is not intended for diagnostic use.

This product is a Human FZD10 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

FZD10

Protein Length

Partial (21-161aa)

Molecular Weight

42.2 kDa

Sequence

21-161aa (Accession # Q9ULW2)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

FZD10

Full Name

Frizzled class receptor 10

Alternative Names

Fz10; FzE7; CD350; FZ-10; hFz10; frizzled-10; frizzled 10, seven transmembrane spanning receptor; frizzled family receptor 10; frizzled homolog 10; FZD10; Frizzled class receptor 10

Gene ID

[11211](#)

UniProt ID

[Q9ULW2](#)