

Product Information

MemDX™ Membrane Protein Human FOLR2 (17-228aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX340**

This product is for research use only and is not intended for diagnostic use.

This product is a Human FOLR2 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

FOLR2

Protein Length

Partial (17-228aa)

Molecular Weight

25.4 kDa

Sequence

17-228aa (Accession # P14207)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

FOLR2

Full Name

Folate receptor beta

Alternative Names

FBP; FOLR1; FR-P3; FRbeta; FR-BETA; BETA-HFR; FBP/PL-1; folate receptor beta; folate receptor 2 (fetal); folate receptor alpha; folate receptor, fetal/placental; folate-binding protein, fetal/placental; placental folate-binding protein

Gene ID

[2350](#)

UniProt ID

[P14207](#)