

Product Information

MemDX™ Membrane Protein Human ALCAM (28-526aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX355**

This product is for research use only and is not intended for diagnostic use.

This product is a Human ALCAM membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

ALCAM

Protein Length

Partial (28-526aa)

Molecular Weight

56.8 kDa

Sequence

28-526aa (Accession # Q13740)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

ALCAM

Full Name

Activated leukocyte cell adhesion molecule

Alternative Names

ALCAM; MEMD; CD166; CD166 antigen; Activated leukocyte cell adhesion molecule

Gene ID

[214](#)

UniProt ID

[Q13740](#)