

Product Information

MemDX™ Membrane Protein Human TNFRSF9 (24-186aa) Expressed in HEK293, C-Mouse Fc and 6×His tag

Cat. No.: **S01YF-0423-KX37**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TNFRSF9 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

TNFRSF9

Protein Length

Partial (24-186aa)

Molecular Weight

54-60 kDa after removal of the signal peptide

Sequence

24-186aa (Accession # Q07011)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Mouse Fc and 6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

TNFRSF9

Full Name

TNF receptor superfamily member 9

Alternative Names

TNFRSF9; ILA; 4-1BB; CD137; CDw137; tumor necrosis factor receptor superfamily member 9; 4-1BB ligand receptor; CD137 antigen; T cell antigen ILA; T-cell antigen 4-1BB homolog; homolog of mouse 4-1BB; induced by lymphocyte activation (ILA); interleukin-activated receptor, homolog of mouse Ly63; receptor protein 4-1BB; TNF receptor superfamily member 9

Gene ID

[3604](#)

UniProt ID

[Q07011](#)