

Product Information

MemDX™ Membrane Protein Human HBEGF (20-160aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX371**

This product is for research use only and is not intended for diagnostic use.

This product is a Human HBEGF membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

HBEGF

Protein Length

Partial (20-160aa)

Molecular Weight

41.9 kDa

Sequence

20-160aa (Accession # Q99075)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

HBEGF

Full Name

Heparin binding EGF like growth factor

Alternative Names

HBEGF; DTR; DTS; DTSF; HEGFL; proheparin-binding EGF-like growth factor; diphtheria toxin receptor (heparin-binding EGF-like growth factor); diphtheria toxin receptor (heparin-binding epidermal growth factor-like growth factor); heparin-binding epidermal growth factor; Heparin binding EGF like growth factor

Gene ID

[1839](#)

UniProt ID

[Q99075](#)