

Product Information

MemDX™ Membrane Protein Human PRLR (25-234aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX423**

This product is for research use only and is not intended for diagnostic use.

This product is a Human PRLR membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

PRLR

Protein Length

Partial (25-234aa)

Molecular Weight

50.5 kDa

Sequence

25-234aa (Accession # P16471)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

PRLR

Full Name

Prolactin receptor

Alternative Names

PRLR; HPRL; MFAB; hPRLrI; RI-PRLR; hPRL receptor; secreted prolactin binding protein; Prolactin receptor

Gene ID

[5618](#)

UniProt ID

[P16471](#)