

Product Information

MemDX™ Membrane Protein Human TNFRSF18 (26-162aa) Expressed in HEK293, C-Mouse Fc and 6×His tag

Cat. No.: **S01YF-0423-KX45**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TNFRSF18 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

TNFRSF18

Protein Length

Partial (26-162aa)

Molecular Weight

50-52 kDa after removal of the signal peptide

Sequence

26-162aa (Accession # Q9Y5U5)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Mouse Fc and 6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

TNFRSF18

Full Name

TNF receptor superfamily member 18

Alternative Names

TNFRSF18; AITR; GITR; CD357; GITR-D; ENERGEN; tumor necrosis factor receptor superfamily member 18; TNF receptor superfamily activation-inducible protein; activation-inducible TNFR family receptor; glucocorticoid-induced TNFR-related protein; TNF receptor superfamily member 18

Gene ID

[8784](#)

UniProt ID

[Q9Y5U5](#)