

Product Information

MemDX™ Membrane Protein Human CDH3 (18-654aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX471**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CDH3 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CDH3

Protein Length

Partial (18-654aa)

Molecular Weight

60.9 kDa

Sequence

18-654aa (Accession # P22223)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

CDH3

Full Name

Cadherin 3

Alternative Names

CDH3; CDHP; HJMD; PCAD; cadherin-3; cadherin 3, type 1, P-cadherin (placental); calcium-dependent adhesion protein, placental; Cadherin 3

Gene ID

[1001](#)

UniProt ID

[P22223](#)