

Product Information

MemDX™ Membrane Protein Human IL10 (19-178aa) Expressed in HEK293, N-Human Fc tag

Cat. No.: **S01YF-0423-KX513**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL10 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL10

Protein Length

Partial (19-178aa)

Molecular Weight

44.8 kDa

Sequence

19-178aa (Accession # P22301)

Product Description

Expression Systems

HEK293

Tag

N-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

IL10

Full Name

Interleukin 10

Alternative Names

CSIF; TGIF; GVHDS; IL-10; IL10A; interleukin-10; T-cell growth inhibitory factor; cytokine synthesis inhibitory factor

Gene ID

[3586](#)

UniProt ID

[P22301](#)