

Product Information

MemDX™ Membrane Protein Human IL15 (49-162aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX516**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL15 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL15

Protein Length

Partial (49-162aa)

Molecular Weight

38.9 kDa

Sequence

49-162aa (Accession # P40933)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

IL15

Full Name

Interleukin 15

Alternative Names

IL-15; interleukin-15

Gene ID

[3600](#)

UniProt ID

[P40933](#)