

## Product Information

### **MemDX™ Membrane Protein Human TLR3 ( ) Expressed in HEK293, C-Human Fc**

Cat. No.: **S01YF-0423-KX553**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TLR3 membrane protein expressed in HEK293.

#### Product Specifications

##### **Host Species**

Human

##### **Target Protein**

TLR3

##### **Molecular Weight**

103.5 kDa

##### **Sequence**

(Accession # O15455)

#### Product Description

##### **Expression Systems**

HEK293

##### **Tag**

C-Human Fc tag

##### **Protein Format**

Soluble

##### **Form**

LYOPH

##### **Purity**

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

##### **Buffer**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

##### **Storage**

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

## Target

### Target Protein

TLR3

### Full Name

Toll like receptor 3

### Alternative Names

TLR3; CD283; IIAE2; Toll like receptor 3

### Gene ID

[7098](#)

### UniProt ID

[Q15455](#)