

Product Information

MemDX™ Membrane Protein Human FLT1 (27-756aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX556**

This product is for research use only and is not intended for diagnostic use.

This product is a Human FLT1 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

FLT1

Protein Length

Partial (27-756aa)

Molecular Weight

108.3 kDa

Sequence

27-756aa (Accession # P17948)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

FLT1

Full Name

Fms related receptor tyrosine kinase 1

Alternative Names

FLT1; FLT; FLT-1; VEGFR1; VEGFR-1; vascular endothelial growth factor receptor 1; fms related tyrosine kinase 1; fms-like tyrosine kinase 1; fms-related tyrosine kinase 1 (vascular endothelial growth factor/vascular permeability factor receptor); tyrosine-protein kinase FRT; tyrosine-protein kinase receptor FLT; vascular permeability factor receptor; Fms related receptor tyrosine kinase 1

Gene ID

[2321](#)

UniProt ID

[P17948](#)