

Product Information

MemDX™ Membrane Protein Human CCR8 (1-35aa) Expressed in HEK293, C-Mouse Fc tag

Cat. No.: **S01YF-0423-KX585**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CCR8 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CCR8

Protein Length

Partial (1-35aa)

Molecular Weight

30.2 kDa

Sequence

1-35aa (Accession # P51685)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Mouse Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CCR8

Full Name

C-C motif chemokine receptor 8

Alternative Names

CY6; TER1; CCR-8; CKRL1; CDw198; CMKBR8; GPRCY6; CMKBRL2; CC-CKR-8; C-C chemokine receptor type 8; CC chemokine receptor 8; CC chemokine receptor CHEMR1; CC-chemokine receptor chemr1; chemokine (C-C motif) receptor 8; chemokine (C-C) receptor 8; chemokine (C-C) receptor-like 2; chemokine receptor-like 1; CCR8; C-C motif chemokine receptor 8

Gene ID

[1237](#)

UniProt ID

[P51685](#)