

Product Information

MemDX™ Membrane Protein Human IL13 (35-146aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX606**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL13 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

IL13

Protein Length

Partial (35-146aa)

Molecular Weight

38.4 kDa

Sequence

35-146aa (Accession # P35225)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

IL13

Full Name

Interleukin 13

Alternative Names

P600; IL-13; interleukin-13

Gene ID

[3596](#)

UniProt ID

[P35225](#)