

Product Information

MemDX™ Membrane Protein Human CDH11 (54-617aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX643**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CDH11 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CDH11

Protein Length

Partial (54-617aa)

Molecular Weight

62.8 kDa

Sequence

54-617aa (Accession # P55287)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

CDH11

Full Name

Cadherin 11

Alternative Names

CDH11; OB; ESWS; CAD11; CDHOB; OSF-4; cadherin-11; cadherin 11, type 2, OB-cadherin (osteoblast); Cadherin 11

Gene ID

[1009](#)

UniProt ID

[P55287](#)