

Product Information

MemDX™ Membrane Protein Human CD70 (63-234aa) Expressed in HEK293, N-Mouse Fc and C-6×His tag

Cat. No.: **S01YF-0423-KX65**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CD70 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CD70

Protein Length

Partial (63-234aa)

Molecular Weight

45.6 kDa

Sequence

63-234aa (Accession # P32970)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

N-Mouse Fc and C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CD70

Full Name

CD70 molecule

Alternative Names

CD70; CD27L; LPFS3; CD27-L; CD27LG; TNFSF7; TNLG8A; CD70 antigen; CD27 ligand; Ki-24 antigen; surface antigen CD70; tumor necrosis factor ligand 8A; tumor necrosis factor ligand superfamily member 7; CD70 molecule

Gene ID

[970](#)

UniProt ID

[P32970](#)