

Product Information

MemDX™ Membrane Protein Human LRP6 (20-1370aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX662**

This product is for research use only and is not intended for diagnostic use.

This product is a Human LRP6 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

LRP6

Protein Length

Partial (20-1370aa)

Molecular Weight

153.1 kDa

Sequence

20-1370aa (Accession # O75581)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

LRP6

Full Name

LDL receptor related protein 6

Alternative Names

ADCAD2; STHAG7; low-density lipoprotein receptor-related protein 6; LRP-6; LRP6; LDL receptor related protein 6

Gene ID

[4040](#)

UniProt ID

[O75581](#)