

Product Information

MemDX™ Membrane Protein Human CLDN6 (29-81aa) Expressed in HEK293, N-Human Fc tag

Cat. No.: **S01YF-0423-KX674**

This product is for research use only and is not intended for diagnostic use.

This product is a Human CLDN6 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

CLDN6

Protein Length

Partial (29-81aa)

Molecular Weight

34.4 kDa

Sequence

29-81aa (Accession # P56747)

Product Description

Expression Systems

HEK293

Tag

N-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CLDN6

Full Name

Claudin 6

Alternative Names

CLDN6; claudin-6; skullin; Claudin 6

Gene ID

[9074](#)

UniProt ID

[P56747](#)