

Product Information

MemDX™ Membrane Protein Human ANGPT2 (275-496aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX679**

This product is for research use only and is not intended for diagnostic use.

This product is a Human ANGPT2 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

ANGPT2

Protein Length

Partial (275-496aa)

Molecular Weight

26.3 kDa

Sequence

275-496aa (Accession # O15123)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

ANGPT2

Full Name

Angiopoietin 2

Alternative Names

ANG2; AGPT2; LMPHM10; angiopoietin-2; Tie2-ligand; angiopoietin-2B; angiopoietin-2a

Gene ID

[285](#)

UniProt ID

[Q15123](#)