

Product Information

MemDX™ Membrane Protein Human LYPD3 (31-286aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX698**

This product is for research use only and is not intended for diagnostic use.

This product is a Human LYPD3 membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

LYPD3

Protein Length

Partial (31-286aa)

Molecular Weight

27.7 kDa

Sequence

31-286aa (Accession # O95274)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

LYPD3

Full Name

LY6/PLAUR domain containing 3

Alternative Names

C4.4A; ly6/PLAUR domain-containing protein 3; 2310061G07Rik; GPI-anchored metastasis-associated protein C4.4A homolog; GPI-anchored metastasis-associated protein homolog; MIG-C4; matrigel-induced gene C4 protein

Gene ID

[27076](#)

UniProt ID

[Q95274](#)