

## Product Information

### **MemDX™ Membrane Protein Human LYPD3 (31-325aa) Expressed in HEK293, C-6×His tag**

Cat. No.: **S01YF-0423-KX699**

This product is for research use only and is not intended for diagnostic use.

This product is a Human LYPD3 membrane protein expressed in HEK293.

#### Product Specifications

##### **Host Species**

Human

##### **Target Protein**

LYPD3

##### **Protein Length**

Partial (31-325aa)

##### **Molecular Weight**

31.7 kDa

##### **Sequence**

31-325aa (Accession # O95274)

#### Product Description

##### **Expression Systems**

HEK293

##### **Tag**

C-6×His tag

##### **Protein Format**

Soluble

##### **Form**

LYOPH

##### **Purity**

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

##### **Buffer**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

##### **Storage**

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

## Target

### Target Protein

LYPD3

### Full Name

LY6/PLAUR domain containing 3

### Alternative Names

C4.4A; ly6/PLAUR domain-containing protein 3; 2310061G07Rik; GPI-anchored metastasis-associated protein C4.4A homolog; GPI-anchored metastasis-associated protein homolog; MIG-C4; matrigel-induced gene C4 protein

### Gene ID

[27076](#)

### UniProt ID

[Q95274](#)