

Product Information

MemDX™ Membrane Protein Human IL15 (49-162aa) Expressed in *E.coli* expression system

Cat. No.: **S01YF-0423-KX727**

This product is for research use only and is not intended for diagnostic use.

This product is a Human IL15 membrane protein expressed in *E.coli* expression system.

Product Specifications

Host Species

Human

Target Protein

IL15

Protein Length

Partial (49-162aa)

Molecular Weight

12.5 KDa

Sequence

49-162aa (Accession # P40933)

Product Description

Expression Systems

E.coli expression system

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Buffer

Lyophilized from a 0.2 µm filtered solution of 20mM PBS, 150mM NaCl, pH 7.0.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

IL15

Full Name

Interleukin 15

Alternative Names

IL-15; interleukin-15

Gene ID

[3600](#)

UniProt ID

[P40933](#)