

Product Information

MemDX™ Membrane Protein Human TNF (35-205aa) Expressed in *E.coli* expression system

Cat. No.: **S01YF-0423-KX802**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TNF membrane protein expressed in *E.coli* expression system.

Product Specifications

Host Species

Human

Target Protein

TNF

Protein Length

Partial (35-205aa)

Molecular Weight

18.8 KDa

Sequence

35-205aa (Accession # P01374)

Product Description

Expression Systems

E.coli expression system

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Buffer

Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

TNF

Full Name

Tumor necrosis factor

Alternative Names

TNF; DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha; tumor necrosis factor; APC1 protein; TNF, macrophage-derived; TNF, monocyte-derived; TNF-a; tumor necrosis factor ligand 1F; tumor necrosis factor ligand superfamily member 2; tumor necrosis factor-alpha; Tumor necrosis factor

Gene ID

[7124](#)

UniProt ID

[P01375](#)