

Product Information

MemDX™ Membrane Protein Cynomolgus ANXA2 (2-339aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX810**

This product is for research use only and is not intended for diagnostic use.

This product is a Cynomolgus ANXA2 membrane protein expressed in HEK293.

Product Specifications

Host Species

Cynomolgus

Target Protein

ANXA2

Protein Length

Partial (2-339aa)

Molecular Weight

39.3 kDa

Sequence

2-339aa (Accession # A0A2K5UQS9)

Product Description

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

ANXA2

Full Name

Annexin A2

Alternative Names

ANXA2; annexin A2

Gene ID

[102139035](#)

UniProt ID

[G7PBL1](#)