

Product Information

MemDX™ Membrane Protein Mouse CD276 (29-244aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX823**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse CD276 membrane protein expressed in HEK293.

Product Specifications

Host Species

Mouse

Target Protein

CD276

Protein Length

Partial (29-244aa)

Molecular Weight

28.7 kDa

Sequence

29-244aa (Accession # Q8VE98)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CD276

Full Name

CD276 antigen

Alternative Names

B7h3; B7RP-2; 6030411F23Rik; CD276 antigen; B7 homolog 3; B7-H3; costimulatory molecule

Gene ID

[102657](#)

UniProt ID

[Q8VE98](#)