

Product Information

MemDX™ Membrane Protein Mouse CTLA4 (36-162aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX825**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse CTLA4 membrane protein expressed in HEK293.

Product Specifications

Host Species

Mouse

Target Protein

CTLA4

Protein Length

Partial (36-162aa)

Molecular Weight

40.1 kDa

Sequence

36-162aa (Accession # P09793)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

CTLA4

Full Name

Cytotoxic T-lymphocyte-associated protein 4

Alternative Names

Cd152; Ly-56; Ctlα-4; cytotoxic T-lymphocyte protein 4; CD152 antigen; cytotoxic T-lymphocyte-associated antigen 4

Gene ID

[12477](#)

UniProt ID

[P09793](#)