

## Product Information

### **MemDX™ Membrane Protein Mouse CD200 (31-232aa) Expressed in HEK293, C-Human Fc tag**

Cat. No.: **S01YF-0423-KX869**

This product is for research use only and is not intended for diagnostic use.

This product is a Mouse CD200 membrane protein expressed in HEK293.

#### Product Specifications

##### **Host Species**

Mouse

##### **Target Protein**

CD200

##### **Protein Length**

Partial (31-232aa)

##### **Molecular Weight**

48.7 kDa

##### **Sequence**

31-232aa (Accession # O54901)

#### Product Description

##### **Expression Systems**

HEK293

##### **Tag**

C-Human Fc tag

##### **Protein Format**

Soluble

##### **Form**

LYOPH

##### **Purity**

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

##### **Buffer**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

**Storage**

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

**Target****Target Protein**

CD200

**Full Name**

CD200 molecule

**Alternative Names**

OX2; Mox2; OX-2 membrane glycoprotein; CD200 antigen; MRC OX-2 antigen; antigen identified by monoclonal antibody MRC OX-2

**Gene ID**

[17470](#)

**UniProt ID**

[Q54901](#)