

Product Information

MemDX™ Membrane Protein Human AFP (19-609aa) Expressed in HEK293, C-Human Fc tag

Cat. No.: **S01YF-0423-KX95**

This product is for research use only and is not intended for diagnostic use.

This product is a Human AFP membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

AFP

Protein Length

Partial (19-609aa)

Molecular Weight

92.6 kDa

Sequence

19-609aa (Accession # P02771)

Product Description

Expression Systems

HEK293

Tag

C-Human Fc tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 90% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target

Target Protein

AFP

Full Name

alpha fetoprotein

Alternative Names

AFPD; FETA; HPAFP; alpha-fetoprotein; alpha-1-fetoprotein; alpha-fetoglobulin

Gene ID

[174](#)

UniProt ID

[P02771](#)