

Product Information

MemDX™ Membrane Protein Human TNF (77-233aa) Expressed in HEK293, C-6×His tag

Cat. No.: **S01YF-0423-KX96**

This product is for research use only and is not intended for diagnostic use.

This product is a Human TNF membrane protein expressed in HEK293.

Product Specifications

Host Species

Human

Target Protein

TNF

Protein Length

Partial (77-233aa)

Molecular Weight

18.1 kDa

Sequence

77-233aa (Accession # P01375)

Product Description

Activity

Yes

Expression Systems

HEK293

Tag

C-6×His tag

Protein Format

Soluble

Form

LYOPH

Purity

Greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Buffer

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.

Storage

The product should be stored at -20°C or lower. Avoid freeze-thaw cycles.

Target**Target Protein**

TNF

Full Name

Tumor necrosis factor

Alternative Names

TNF; DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha; tumor necrosis factor; APC1 protein; TNF, macrophage-derived; TNF, monocyte-derived; TNF-a; tumor necrosis factor ligand 1F; tumor necrosis factor ligand superfamily member 2; tumor necrosis factor-alpha; Tumor necrosis factor

Gene ID

[7124](#)

UniProt ID

[P01375](#)